

# Samsung PT10V

## Profilübersicht



Parameter	Comprehensive Profil	Electrolyte Profil	Kidney Profil	Large Animal Profil	Liver Profil	Pre-Surgical Profil
A/G ratio*	X			X	X	X
ALB	X		X	X	X	X
ALP	X			X	X	X
ALT	X				X	X
AMY	X		X			
AST				X	X	
B/C ratio*	X		X	X		X
BUN	X		X	X	X	X
CA	X		X			
CHOL	X					
CK				X		
CL		X				
CREA	X		X	X		X
GGT				X	X	
GLOB*	X			X	X	X
GLU	X		X	X	X	X
K		X				
NA		X				
NA/K ratio*		X				
PHOS	X		X			
TBIL	X			X	X	
TP	X			X	X	X
TRIG	X					

\* berechnete Parameter

scil animal care company GmbH  
 Dina-Weissmann-Allee 6  
 68519 Viernheim  
 Tel.: +49 (0) 6204 78 90 222  
 Fax: +49 (0) 6204 78 90 200  
 Email: info-de@scilvet.com  
 www.scilvet.de

